

**DNase I-RNase Free**

SB-ME-002S- 10,000 U (10 U/ul)

SB-ME-002L-50,000 U (10 U/ul)

**Store at -20°C****Description**

This product is a recombinant form of DNase I from bovine pancreas, expressed in *Pichia pastoris*. DNase I is a DNA-specific endonuclease that hydrolyses the phosphodiester linkages of double-stranded or single-stranded DNA to a mixture of oligo- and mononucleotides. The enzyme requires divalent cations for maximal activity. It is a glycoprotein of a molecular weight of ~39KD. The recombinant enzyme is produced without using any animal cells or other materials derived from animals.

Characteristics	Applications	Quality control
<ul style="list-style-type: none"> <li>- Molecular weight : ~39 Kda</li> <li>- DNase I is heterogeneously N-glycosylated, so it appears as two bands in gel electrophoresis.</li> </ul>	<ul style="list-style-type: none"> <li>- Remove genomic DNA from RNA preparations prior to RT-PCR</li> <li>- Isolate DNA-free RNA after in vitro transcription reactions</li> <li>- Map DNase-sensitive regions in eukaryotic DNA</li> <li>- Nick translation</li> </ul>	<ul style="list-style-type: none"> <li>- Purity : &gt;95% on SDS-PAGE</li> <li>- RNase-free</li> <li>- Protease-free</li> </ul>

**Storage buffer**20 mM Tris-HCl, 50 mM NaCl, 2 mM CaCl<sub>2</sub> , 2 mM MgCl<sub>2</sub> ,1 mM DTT, 0.1 mg/ml Pefabloc SC, 50% Glycerol, pH 7.6**Diluent Buffer**

25 mM Tris-HCl, 50% Glycerol, pH 7.6

**Heat inactivation**

Inactivated by heating at 75°C for 10 min in the presence of EDTA (final 8 mM)

Note DNase I requires divalent cations for maximal activity. The enzyme is inhibited by metal chelating agents like EDTA. DNase I can be inactivated and removed by phenol extraction according to standard protocols.

**Limitation of Liability**

In no event will SARD BIOSCIENCES be liable for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

**Standard reaction conditions**

<b>DNA</b>	<b>2µg</b>
10X DNase I buffer	2 µl
DNase I, RNase-free	1~2 U
Water, RNase-free	up to 20 µl

→ Incubate at 37°C for 10 min

'Standard conditions' is only a general recommendation. The experimental conditions should be adjusted according to the purpose and sample